Angiotensin Converting Enzyme (ACE)

CPT Code: 82164 Order Code: 683

ABN Requirement: No

Alternative Name(s): Angiotensin-1-Converting Enzyme, ACE

Specimen:

Preferred: Serum

Alternate: Lithium Heparin Plasma

Volume: 1.0 mL

Minimum Volume: 0.5 mL

Container:

Preferred: Serum Separator Tube (SST; Gel Barrier) Plasma collected in: Sodium

Heparin (green-top) tube Lithium Heparin (green-top) tube

Alternate: Green Top (Lithium Heparin) Tube

Collection:

Serum (preferred):

- 1. Collect and label sample according to standard protocols.
- 2. Gently invert tube 5 times immediately after draw. Do not shake.
- 3. Allow blood to clot 30 minutes.
- 4. Centrifuge for 10 minutes.

Lithium Heparin Plasma (alternative):

- 1. Collect and label sample according to standard protocols.
- 2. Gently invert tube 10 times immediately after draw. Do not shake.
- 3. Centrifuge for 10 minutes.
- 4. Aliquot plasma into transport tube labeled as "Lithium Heparin Plasma" and cap tightly.

Discard original tube.

Collection Instructions: Centrifuge specimen within 1 hour of collection.

Transport: Store specimen at 2°C to 8°C after collection and ship the same day

per packaging instructions included with the provided shipping box.

Stability:

Ambient (15-25°C): 7 days **Refrigerated (2-8°C):** 7 days

Frozen (-20°C): 60 days

Causes for Rejection: Gross hemolysis; Grossly lipemic

Methodology: Spectrophotometry (SP)

Turn Around Time: 2 to 4 days

Reference Range:

Age	Reference Range (U/L)
6 Months-17 Years	13-100 U/L
Adult	9-67 U/L

Clinical Significance: This test is useful in evaluating patients presenting with sarcoidosis, Gaucher's disease, and lymphangiomyomatosis in that the enzyme is increased in these clinical settings.

The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.