## Lipoprotein Fractionation, NMR

**CPT Code:** 83704 **Order Code:** 37847

Includes: LDL Particle Number (LDL-P), Small LDL Particle Number (Small LDL-P), LDL Size, HDL Particle Number (HDL-P), Large HDL Particle Number (Large

HDL-P), HDL Size, Large VLDL Particle Number (Large VLDL-P), VLDL Size

**ABN Requirement:** No

**Specimen**: Serum **Volume**: 2.0 mL

Minimum Volume: 1.0 mL

Container: Red Top (no gel barrier) tube

## **Collection**:

1. Collect and label sample according to standard protocols.

- 2. Gently invert tube 5 times immediately after draw. DO NOT SHAKE.
- 3. Allow blood to clot 30-60 minutes.
- 4. Centrifuge for 10 minutes.
- 5. Aliquot serum into a labeled transport tube labeled as "Red Top Serum" and cap tightly.
- 6. Store and transport refrigerated. It is acceptable to store frozen and ship refrigerated. Do not ship frozen on dry ice.

Fasting: Patient should be fasting 12 hours.

**Transport:** Store serum at 2°C to 8°C after collection and ship the same day per packaging instructions included with the provided shipping box.

## **Stability:**

Ambient (15-25°C): Unacceptable

Refrigerated (2-8°C): 7 days

**Frozen (-20°C):** 10 days

**Deep Frozen (-70°C):** 28 days

Causes for Rejection: Specimens other than serum; improper labeling; samples not stored properly; samples older than stability limits; gross hemolysis; gross

lipemia; serum collected using a gel barrier tube

Methodology: Nuclear Magnetic Resonance (NMR)

**Turn Around Time:** 3 days

## **Relative Risk Range:**

**Clinical Significance:** The Lipoprotein Fractionation NMR test is used to help assess the risk for cardiovascular disease (CVD) in patients with intermediate or high risk based on traditional or emerging risk factors, and to assess therapeutic response in patients undergoing lipid-lowering therapy, by quantification of the number and size of lipoprotein particles.

The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.