

# Zinc, RBC

**CPT Code:** 84630

**Order Code:** 6354

**ABN Requirement:** No

**Specimen:** Red blood cells; see collection instructions

**Volume:** 0.5 mL

**Minimum Volume:** 0.3 mL

**Container:**

**Preferred:** EDTA trace metal-free (Royal Blue top) Tube

**Alternative:** EDTA (Lavender Top) Tube, Sodium Heparin (Royal Blue Top) Tube, Sodium Heparin lead-free (Tan Top) Tube, Sodium or Lithium Heparin (Green Top) Tube

**Collection:**

**Red blood cells collected in trace-metal free tube:** Centrifuge to separate plasma from RBCs and discard plasma. Use the EDTA royal blue-top trace metal evacuated tube for RBC testing.

**Red blood cells collected in other tubes:** Centrifuge to separate plasma from RBCs and discard plasma. Transfer RBCs to a plastic transfer vial, from a Quest Diagnostics trace element and metal free collection kit.

**Collection Instructions:** Carefully clean skin prior to venipuncture. Avoid hemolysis. Avoid worksite collection.

**Fasting:** Fasting is not required.

**Transport:** Store serum at 2°C to 8°C after collection and ship the same day per packaging instructions included with the provided shipping box.

**Stability:**

**Ambient (15-25°C):** 7 days

**Refrigerated (2-8°C):** 10 days

**Frozen (-20°C):** Not Acceptable

**Causes for Rejection:** Improper labeling; samples not stored properly; samples

older than stability limits; clotted; moderate hemolysis

**Methodology:** Inductively Coupled Plasma/Mass Spectrometry (ICP/MS)

**Turn Around Time:** 4 to 6 days

**Reference Range:**

	mg/L
All Ages	9.0-14.7

**Clinical Significance:** Zinc is an essential element involved in a myriad of enzyme systems including wound healing, immune function, and fetal development. Zinc measurements are used to detect and monitor industrial, dietary, and accidental exposure to zinc. Also, Zinc measurements may be used to evaluate health and monitor response to treatment.

*The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.*